


The **Peripheral** Edge in Every **Flex**




MicroFlex

Cobalt Chromium Peripheral Stent System




COBALT CHROMIUM STENT WITH 'C' & 'S' SHAPED POLYLINK

- CoCr alloy provide superior radial strength and excellent radiopacity.
- 14 crown stent design for uniform and optimum vessel scaffolding.



14-22 atm RBP DELIVERS HIGH PRESSURE CAPABILITY ON LOW PROFILE SYSTEM



PRECISION DEPLOYMENT WITH MINIMAL RECOIL & FORESHORTENING

- › Designed for precise expansion with minimal recoil and foreshortening.
- › Predictable performance supports precise stent placement across various lesion sites.

ULTRA-LOW PROFILE DESIGN FOR OPTIMIZED TRACKABILITY

- › Low profile shaft facilitates smooth advancement through tortuous anatomies.
- › Excellent compatibility with standard delivery systems, reducing procedural complexity.

INDICATIONS



Peripheral
Artery Diseases



Renal Artery



Palliation in Malignant
Stricture of the Biliary Tree.

SPECIFICATION

Stent Material	CoCr
Stent Diameter	5, 6, 7, 8, 9, 10
Stent Length	17, 27, 37, 47, 57, 67
Strut Thickness	0.120 mm
System Type	OTW
Nominal Pressure	8-10 atm
Rated Burst Pressure	14-22 atm
Shaft Outer Diameter	5 Fr
Radiopaque Marker	2-Marker of Pt/Ir
Guidewire Compatibility	0.035"
Sheath Compatibility	6 / 7 Fr

ORDERING INFORMATION

Usable Catheter Length 135 cm (OTW 0.035")						
Stent Length (mm)						
Diameter (mm)	17	27	37	47	57	67
5	MF35A0517	MF35A0527	MF35A0537	MF35A0547	MF35A0557	MF35A0567
6	MF35A0617	MF35A0627	MF35A0637	MF35A0647	MF35A0657	MF35A0667
7	MF35A0717	MF35A0727	MF35A0737	MF35A0747	MF35A0757	MF35A0767
8	MF35A0817	MF35A0827	MF35A0837	MF35A0847	MF35A0857	MF35A0867
9	MF35A0917	MF35A0927	MF35A0937	MF35A0947	MF35A0957	MF35A0967
10	MF35A1017	MF35A1027	MF35A1037	MF35A1047	MF35A1057	MF35A1067

Corporate Office - Sahajanand Laser Technology Ltd.

E-30, Electronic Estate, G.I.D.C., Sector-26, Gandhinagar - 382 028, Gujarat, India.

Tel : +91-79-6170-1616 | +91-99-2503-6495

Email : mkt@SLTL.com

www.SLTLmedical.com

Manufacturing Site - Sahajanand Laser Technology Ltd.

A-8, Electronic Estate, G.I.D.C., Sector - 25, Gandhinagar - 382 016, Gujarat, India.

Tel : +91-99-2503-6495

MF/CAT/00



Caution: Indication, Contraindication, Warnings, Precautions and other instructions for use can be found in the product labeling supplied with each device.