


zephyr[®]_M

Cobalt Chromium Peripheral Stent System

Featuring an innovative overexpansion stent design,
it ensures optimal vessel fit and reliable post-procedural expansion performance.

Features:

- Unique 12-crown design provides optimal radial strength and scaffolding.
- “S” and “C” shaped polylinks enhance flexibility and conformability.
- Minimizes foreshortening for accurate stent placement.



Technical Specifications:

Stent Material CoCr	Strut Thickness 0.120 mm	Delivery System OTW	Guidewire Compatibility 0.035"
Foreshortening <1 %	Stent Recoil <6 %	Usable Device Length 135 cm	Min. Recommended Sheath 6F

Ordering Information:

Clinical Indications:

- Treatment of de novo or restenotic atherosclerotic lesions in peripheral arteries
- Treatment of native and /or recurrent coarctation of aorta
- Stenosis of main / branch pulmonary artery

Length	Diameter					
	5	6	7	8	9	10
13	ZPMB0513	ZPMB0613	ZPMB0713	ZPMB0813	ZPMB0913	ZPMB1013
16	ZPMB0516	ZPMB0616	ZPMB0716	ZPMB0816	ZPMB0916	ZPMB1016
20	ZPMB0520	ZPMB0620	ZPMB0720	ZPMB0820	ZPMB0920	ZPMB1020
23	ZPMB0523	ZPMB0623	ZPMB0723	ZPMB0823	ZPMB0923	ZPMB1023
25	ZPMB0525	ZPMB0625	ZPMB0725	ZPMB0825	ZPMB0925	ZPMB1025
30	ZPMB0530	ZPMB0630	ZPMB0730	ZPMB0830	ZPMB0930	ZPMB1030
34	ZPMB0534	ZPMB0634	ZPMB0734	ZPMB0834	ZPMB0934	ZPMB1034
37	ZPMB0537	ZPMB0637	ZPMB0737	ZPMB0837	ZPMB0937	ZPMB1037
40	ZPMB0540	ZPMB0640	ZPMB0740	ZPMB0840	ZPMB0940	ZPMB1040
43	ZPMB0543	ZPMB0643	ZPMB0743	ZPMB0843	ZPMB0943	ZPMB1043
47	ZPMB0547	ZPMB0647	ZPMB0747	ZPMB0847	ZPMB0947	ZPMB1047
50	ZPMB0550	ZPMB0650	ZPMB0750	ZPMB0850	ZPMB0950	ZPMB1050
53	ZPMB0553	ZPMB0653	ZPMB0753	ZPMB0853	ZPMB0953	ZPMB1053
57	ZPMB0557	ZPMB0657	ZPMB0757	ZPMB0857	ZPMB0957	ZPMB1057
60	ZPMB0560	ZPMB0660	ZPMB0760	ZPMB0860	ZPMB0960	ZPMB1060
65	ZPMB0565	ZPMB0665	ZPMB0765	ZPMB0865	ZPMB0965	ZPMB1065
Compliance Chart	NP:10 atm	NP:10 atm	NP:10 atm	NP:10 atm	NP:8 atm	NP:8 atm
	RBP: 22 atm	RBP: 22 atm	RBP: 20 atm	RBP: 18 atm	RBP: 18 atm	RBP: 14 atm

Corporate Office - Sahajanand Laser Technology Ltd.
E-30, Electronic Estate, G.I.D.C., Sector-26,
Gandhinagar - 382 028, Gujarat, India.
Tel : +91-79-6170-1616
Cell : +91-99-250-36495
Email : mkt@SLTL.com

Manufacturing Site - Sahajanand Laser Technology Ltd.
A-8, Electronic Estate, G.I.D.C., Sector - 25,
Gandhinagar - 382 016, Gujarat, India.
Cell : +91-99-250-36495

ZM/CAT/00

www.SLTLmedical.com
www.SLTL.com