


Understanding the Art of Stenting



Easyflex

EasyFlex Stainless Steel Coronary Stent System




One of the
most
flexible stent
designs

Easyflex

EasyFlex Stainless Steel Coronary Stent System

Two opposite articulations between every stent segment guarantee multidirectional flexibility



- Easy side branch access
- Safe and accurate delivery
- 316 LVM stainless steel
- Minimal foreshortening
- Unique star segments
- Low crimping profile
- Biocompatible
- Balloon expandable
- High flexibility
- Good trackability

Stent Material
316LVM Stainless Steel

Strut Thickness
90 μ

Mean Recoil
2.30%

Mean Foreshortening
0.60%

Crossing Profile
0.97 mm (3mm)

Balloon Overhang
<0.5mm

Usable Catheter Length
145cm

Guidewire Compatibility
0.014"

Guide catheter Compatibility
5F

MRI Compatibility
MR Conditional

Compliance
Nominal Pressure: 6 bar
RBP: 14-16 bar

Easyflex

EasyFlex Stainless Steel Coronary Stent System

CATALOG NO :

Length	Diameter (Ø)					
	2.25mm	2.50mm	2.75mm	3.00mm	3.50mm	4.00mm
07mm	EF 0722	EF 0725	EF 0727	EF 0730	EF 0735	EF 0740
10mm	EF 1022	EF 1025	EF 1027	EF 1030	EF 1035	EF 1040
13mm	EF 1322	EF 1325	EF 1327	EF 1330	EF 1335	EF 1340
15mm	EF 1522	EF 1525	EF 1527	EF 1530	EF 1535	EF 1540
17mm	EF 1722	EF 1725	EF 1727	EF 1730	EF 1735	EF 1740
20mm	EF 2022	EF 2025	EF 2027	EF 2030	EF 2035	EF 2040
24mm	EF 2422	EF 2425	EF 2427	EF 2430	EF 2435	EF 2440
28mm	EF 2822	EF 2825	EF 2827	EF 2830	EF 2835	EF 2840
33mm	EF 3322	EF 3325	EF 3327	EF 3330	EF 3335	EF 3340
38mm	EF 3822	EF 3825	EF 3827	EF 3830	EF 3835	EF 3840
42mm	EF 4222	EF 4225	EF 4227	EF 4230	EF 4235	EF 4240
45mm	EF 4522	EF 4525	EF 4527	EF 4530	EF 4535	EF 4540

International Head Quarter - Sahajanand Laser Technology Ltd.

E-30, G.I.D.C., Electronic Estate, Sector-26,
Gandhinagar - 382 028. (Gujarat) India.

Tel. : +91 - 79 - 2324 1379 / 80,

Fax : +91 - 79 - 2324 1465

Email : mkt@SLTL.com

Log on to :

www.SLTLMedical.com

www.SLTL.com

India Marketing - Sahajanand Laser Technology Ltd.

41, Newyork Tower-A, Nr. Thaltej Cross Road,
S.G. Highway, Ahmedabad - 380 054. Gujarat, India.

Tel.: +91-79-2685 4537 / 38 , Fax: +91-79-2685 4539

Manufacturing - Sahajanand Laser Technology Ltd.

A-8, G.I.D.C., Electronic Estate, Sector-25,
Gandhinagar - 382 016. (Gujarat) India.

Tel.: +91 - 79 - 2328 7462, Fax: +91 - 79 - 2324 1465



CAUTION: Indication, Contraindication, Warnings, Precautions and other instructions for use can be found in the product labeling supplied with each device.
EasyFlex is a registered trademark of Sahajanand Laser Technology Ltd. © 2020 SLTL. All right reserved.