


**THE ULTIMATE CHOICE FOR  
PERIPHERAL ARCHITECTURE**



**zephyr<sub>L</sub>**

**Cobalt Chromium Peripheral Stent**

## CoCr STENT DESIGNED FOR LARGE VESSEL SUPPORT

- US patented unique zigzag hybrid stent design combined with “S” & “C” shaped polylink
- Cobalt chromium material imparts low crimping profile, high radial strength and increased radiopacity



## SPECIFICATIONS:

Stent Material	Strut Thickness	Surface / Artery Ratio	Recommended Sheath	Foreshortening
<b>Cobalt Chromium</b>	<b>0.009"</b>	<b>18%</b>	<b>8 F (+ 1F)</b> (depending on the PTA)	<b>&lt; 5%</b>
Stent Design <b>Zigzag Hybrid</b>	Links <b>Flexible Polylinks</b>	Stent Recoil <b>&lt; 3%</b>	Dilatation Range <b>Upto 18 mm</b>	Stent Diameter Unexpanded <b>4.7 mm</b>

## ORDERING INFORMATION:

	Length (mm)											
	14	18	23	28	33	38	43	48	53	58	63	68
<b>Product Codes</b>	ZPLG14	ZPLG18	ZPLG23	ZPLG28	ZPLG33	ZPLG38	ZPLG43	ZPLG48	ZPLG53	ZPLG58	ZPLG63	ZPLG68

\* SLTL Recommends use of BIB Balloon or Cristal Balloon for Zephyr stent family

**Corporate Office - Sahajanand Laser Technology Ltd.**  
 E-30, Electronic Estate, G.I.D.C., Sector-26,  
 Gandhinagar - 382 028, Gujarat, India.  
 Tel : +91-79-2328-7461-63-64, Fax : +91-79-2328-7470  
 Cell : +91-90-9904-0080  
 Email : mkt@SLTL.com  
**Log on to :**  
[www.SLTLMedical.com](http://www.SLTLMedical.com)  
[www.SLTL.com](http://www.SLTL.com)

**Manufacturing Site - Sahajanand Laser Technology Ltd.**  
 A-8, Electronic Estate, G.I.D.C., Sector-25,  
 Gandhinagar - 382 016, Gujarat, India.  
 Cell : +91-90-9904-0090, Fax : +91-79-2328-7470